

PRODUCT SPECIFICATION

SecondSkin

POWERED BY
SEMPERIT

Light

Spec code: N-1135-24B
Product name: SecondSkin Light
Date of issue: 28.06.2016

GENERAL INFORMATION

Type: single-use disposable examination or protective glove
Marking: all information on dispenser box
Shape: ambidextrous - straight fingers
Material: Nitrile Butadiene Rubber (NBR)
Colour: lavender-blue
Inner: online chlorinated / powderfree
Outer: no treatment
Cuff/Surface: rolled cuff / finger textured

Thickness: <small>(double measured)</small>	Finger: min. 0.14 mm	Palm: min. 0.12 mm	Cuff: min. 0.10 mm
---	-----------------------------	---------------------------	---------------------------


Residual powder: ≤ 2 mg / glove	Force at break <small>(according to EN 455-2)</small>	med. 6.0 N <small>(during shelf life)</small>
--	---	--

Length:	median 240 mm
----------------	---------------

Available Sizes:	S (6-7) M (7-8) L (8-9) XL (9-10)
-------------------------	-----------------------------------

Storage: 3 years storage life (date of expiration see on the box). Keep area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolour the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead the accelerated aging and should be avoided.

STANDARDS, GUIDELINES AND QUALITY CERTIFICATION




Quality certification: ISO 9001, ISO 13485
Conformity to guidelines: Personal Protective Equipment 89/686/EEC: CE cat. III
 Medical Device Directive 93/42/EEC: CE Class I
Conformity to standards: EN 420 EN 374-2: specified EN 374-3: Class 6
 EN 455 level 2 - AQL 1.5 Sodium Hydroxide 40%
 ASTM D6319 level 3: AQL <0.65 (G1)
 level 2: AQL <1.5 (G1)
 level 1: <4.0 (S4)
Conformity to food regulation: 1935/2004 

WARNING AND REPORTING SYSTEM

Potentially allergic ingredients: Dihiocarbamate type, Zinc-mercaptobenzothiazol
(For further information, the actual version of the ingredients list is given on request.)

Medical Device Villigance and Reporting System: According to the official reporting criteria of MDD, incidents caused by examination gloves or surgical gloves must be reported immediately to our medical device reporting officer. Tel. +43 2630 310 ext. -510; Fax: 549

Caution! This product contains accelerators which may cause allergic reactions!

Established (Product Manager):  N. Vuksanovic	Reviewed (Sales Director):  G. Cappra	Approved (CEO):  N. Saladin
---	---	---

Edi Baur AG
 Im Schrännbrunnen 1, 8320 Fehraltorf, Switzerland
 Tel.: +41 43 355 70 20 / Fax: +41 43 355 70 30
 info@edibaur.com - www.edibaur.com

SecondSkin Light
 N-1135-24B
 28.06.2016